

**Product: Expression Arrest™ pSM2 empty vector**

Catalog #: RHS1763


The laboratory of Dr. Greg Hannon at Cold Spring Harbor Laboratory (CSHL) has created an RNAi Clone Library comprised of multiple short-hairpin RNAs (shRNAs) specifically targeting annotated human and mouse genes. The shRNA Library permits rapid, cost efficient, loss-of-function genetic screens and rapid tests for genetic interactions to be performed in mammalian cells. Each shRNA has been cloned and sequence verified to ensure a match to the target gene.

The shRNA Control Vector (See Figure 1) is a negative control for any transfection experiment performed using the Expression Arrest™ shRNA-containing expression vectors. This control vector is pSHAG-MAGIC 2 without an insert.

Each vial of the shRNA Control Vector is shipped at a concentration of 0.25µg/µl in a total volume of 40µl, thus providing a total amount of 10µg vector DNA.

**shRNA vector storage**

The vector DNA is shipped in a microfuge tube at room temperature and should be stored at -20°C or -80°C.



**Figure 1: Vector Map of pSHAG-MAGIC 2**



*Note: The pSM vectors must be transformed into PIR1 competent bacteria. These plasmids harbor a conditional bacterial origin of replication, which requires the expression of the “pir1” gene to be rendered functional.*

### **Antibiotic Resistance**

pSHAG-MAGIC 2 contains 3 antibiotic resistance markers. (See Table 1)

Table 1: Antibiotic Resistance Conveyed by pSHAG-MAGIC 2

<b>Antibiotic</b>	<b>Concentration</b>	<b>Utility</b>
Chloramphenicol	25µg/ml	Bacterial selection marker (shRNA insert)
Kanamycin	25µg/ml	Bacterial selection marker (vector)
Puromycin		Mammalian selectable marker

The pSHAG-MAGIC expression vector is a self-inactivating (SIN) MSCV retroviral vector, containing a NheI/XbaI deletion in the U3 region of the 3' LTR. Self-inactivating MSCV particles can be produced by transfection into commonly available retroviral packaging lines.

### **Useful References:**

Carmell, M.A., Zhang, L., Conklin, D.S., Hannon, G.J., Rosenquist, T.A., *Germline transmission of RNAi in mice*. Nat Struct Bio 10, 91-92. (2003).

Conklin, D., Hannon, G., Bernstein, E., Caudy, A., Paddison, P., *Short Hairpin RNAs (shRNAs) Induce Sequence-Specific Silencing in Mammalian Cells*. Genes & Development. 16:948-958. (2002).

Hannon, G., Caudy, A., Paddison, P., *Stable Suppression of Gene Expression by RNAi in Mammalian Cells*. PNAS 99(3):1443-1448 (2003).

Hemann, M.T., Fridman, J.S., Zilfou, J.T., Hernando, E., Paddison, P., Cordon-Cardo, C., Hannon, G.J., Lowe, S.W. *An epi-allelic series of p53 hypomorphs created by stable RNAi produces distinct tumor phenotypes in vivo*. Nat Genet 33, 396-400. (2003).

McCaffrey A.P., Meuse L., Pham T.T., Conklin D.S., Hannon G.J., Kay M.A. *RNA interference in adult mice*. Nature (418) 38-39. (2002).

Paddison, P.J., Silva, J.M., Conklin, D.S., Schlabach, M., Li, M., Aruleba, S., Balija, V., O'Shaughnessy, A., Gnoj, L., Scobie, K., Chang, K., Westbrook, T., Sachidanandam, R., McCombie W.R., Elledge, S.J., & Hannon, G.J. *A resource for large-scale RNAi based screens in mammals*. Nature, March 25; 427-431. (2004)



Paddison, P.J., Caudy, A.A., Bernstein, E., Hannon, G.J. & Conklin, D.S. *Short hairpin RNAs (shRNAs) induce sequence-specific silencing in mammalian cells.* Genes Dev 16, 948-58. (2002).

Taira, K., Miyagishi, M. *U6 promoter-driven siRNAs with four uridine 3' overhangs efficiently suppress targeted gene expression in mammalian cells.* Nature Biotechnology 19, 497-500. (2002).

**Limited Use License:**

*This product is covered by several patent applications owned by Cold Spring Harbor Laboratory. The purchase of this product conveys to the buyer the limited, non-exclusive, non-transferable right (without the right to resell, repackage, or further sublicense) under these patent rights to perform the RNAi knockdown methods using the RNAi inducing vectors claimed in those patent applications for research purposes solely in conjunction with this product. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of this product does not include nor carry any right or license to use, develop, or otherwise exploit this product commercially, and no rights are conveyed to the buyer to use the product or components of the product for any other purposes, including without limitation, provision of services to a third party, generation of commercial databases, or clinical diagnostics or therapeutics.*

*In addition, any commercial organization that purchases or desires to purchase RNAi clones to greater than one thousand (1000) different genes may be outside the above research license and will contact Cold Spring Harbor for a license.*

*This product is sold pursuant to a license from CSHL, and CSHL reserves all other rights under these patent rights. For information on purchasing a license to the patent rights for uses other than in conjunction with this product or to use this product for purposes other than research, please contact the CSHL Office of Technology Transfer at, 516-367-8312.*